



Infuse, disconnect, & discard.

1. Ensure clamps are open (Homepump Eclipse* device tubing as well as on patient's access device).
2. Check to make sure tubing is free of kinks.
3. Ensure the Homepump Eclipse* pump operates at the correct temperature:
 - a. Pump must be at room temperature before infusion.
 - b. Pump and tubing must be outside patient's clothing.

How much time does it take for the pump to reach room temperature?

Size (ml)	Remove from Refrigerator (hours)	Remove from Freezer (hours)
50	7	12
100	8	16
250	12	25
400	15	29
500	18	30

Tip: Advise patient to take the pump with the next dose of medicine out of the refrigerator while a dose is infusing. This will allow each dose enough time to reach room temperature prior to use.

How do I know the pump is infusing?

Depending on the size of the pump, a change in the size and appearance of the pump should be apparent fairly quickly (30 minutes or less).

When is the infusion is complete?

Infusion times vary based on the size and flow rate of the pump. When the delivery time has passed and the pump is no longer inflated, the infusion is complete. A hard tube will be present in the middle of the pump.

There are inherent risks in all medical devices. Please refer to the product labeling for **Indications, Cautions, Warnings and Contraindications**. Failure to follow the product labeling could directly impact patient safety. Physician is responsible for prescribing and administering medications per instructions provided by the drug manufacturer. Refer to www.avanosmedicaldevices.com for product safety Technical Bulletins.

Rx only. *Registered Trademark or Trademark of Avanos Medical, Inc. or its affiliates. © 2015 AVNS. All rights reserved.
MK-00298_Rev2 03/2015