

MORE SIZES PROVIDE BETTER OPTIONS

TAKE YOUR PATIENT CARE TO
GREATER LENGTHS

21
MORE SIZES
THAN AMT®¹

**NOW
UP TO
7.0cm**
STOMA
LENGTHS²


THE LARGEST
low-profile GJ
offering available.¹

**LOW-PROFILE HEAD
FOR PATIENT COMFORT**
Separate jejunal and gastric
ports for simultaneous feeding
and gastric decompression³


**TUNGSTEN WEIGHTED
IN ALL SIZES**
For increased radiopacity
and designed to help
maintain positioning⁴

**ONLY 1 TYPE OF
EXTENSION
SET REQUIRED⁶**

**MULTIPLE LARGE
FEEDING EXIT PORTS**
Designed to improve flow
and minimize clogging⁵



**SMALLER FEED HEAD FOOTPRINT⁷
THAN AMT® G-JET**





THE COMPLETE PORTFOLIO

From your youngest patient to your oldest patient, *Avanos has you covered.*
 We have the MOST comprehensive low-profile feeding tube portfolio on the market to meet your clinical needs⁸.

		Jejunal Length (cm)			
		15	22	30	45
Stoma Length (cm)	1.0	14 Fr 16 Fr	14 Fr 16 Fr	N/A	N/A
	1.2	14 Fr 16 Fr	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr	N/A
	1.5	14 Fr 16 Fr	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr 22 Fr
	1.7	14 Fr 16 Fr	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr 22 Fr
	2.0	N/A	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr 22 Fr
	2.3	N/A	N/A	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr 22 Fr
	2.5	N/A	N/A	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr 22 Fr
	2.7	N/A	N/A	14 Fr 16 Fr 18 Fr	14 Fr 16 Fr 18 Fr 22 Fr
	3.0	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
	3.5	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
	4.0	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
	4.5	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
	5.0	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
	5.5	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
	6.0	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
	6.5	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr
7.0	N/A	N/A	N/A	14 Fr 16 Fr 18 Fr 22 Fr	

How To Configure Your MIC-KEY* GJ Feeding Tube Kit Product Code:

PLEASE NOTE: Codes beginning in 0270 are Non-ENFit® codes and codes beginning in 8270 are ENFit® codes.

NON-ENFit® CODE EXAMPLE: 0270 + _____ (Fr) + _____ (cm) + _____ (cm) = Product Code: _____
French Size Stoma Length Jejunal Tail Length NON-ENFit® CODE Example: 0270-14-1.0-15
(Outer Diameter)

ENFit® CODE EXAMPLE: 8270 + _____ (Fr) + _____ (cm) + _____ (cm) = Product Code: _____
French Size Stoma Length Jejunal Tail Length ENFit® CODE Example: 8270-14-1.0-15
(Outer Diameter)

1. Data on file (CLAIM 07722). 2. Data on file (CLAIM 07723). 3. Data on file (CLAIM 07429). 4. Data on file (CLAIM 07725).
 5. Data on file (CLAIM 07724). 6. Data on file (CLAIM 07356). 7. Data on file (CLAIM 07355). 8. Data on file (CLAIM 07748)

For more information please visit: avanosmedicaldevices.com Call 1-844-4AVANOS (1-844-428-2667) in the United States and Canada.