




Medication Filling & Administration


from Hospital to Home

NEOMED[®]
AN **AVANOS** COMPANY


Medication Filling & Administration




Syringes
Page 4




DoseMate & DoseMate DL
Page 8



Blunt Needle
Page 10




Tamper-Evident Caps
Page 6




Syringe-to-Syringe Coupler
Page 9




Fill Cap Coupler
Page 10



Pharmacy Caps
Page 7



Transfer Straws
Page 9



NeoConnect at Home
Page 11

Join the global patient safety movement to **ENFit**

ENFit is a unique enteral-specific design that mandates a specific direction of fluid flow (female to male) by providing a secure, locking connection to minimize inadvertent disconnections. ENFit is designed to support patient safety initiatives and mitigate the risk of misconnection with non-enteral devices.

How have **syringes** changed?



Legacy



ENFit

Syringes to administer medicine, flush, hydrate, or bolus feed through enteral tubes now require a precise enteral-specific fitment. To ensure small volume dosing accuracy, GEDSA requires an ENFit Low Dose Tip for syringe sizes of 5 mL or smaller.



For accurate enteral dosing of small doses, syringe sizes of 5 mL or smaller require an ENFit Low Dose Tip (LDT) syringe design.




“Adapters for syringes add risk”

“The add-on devices appear to have **significant dead space** that will allow the accumulation of fluid during administration, which will never reach the patient. Thus, these **adapters can cause inaccurate liquid dosing of small volume liquids** in neonates and pediatric patients, and in adults who are receiving drugs that have a narrow therapeutic index.”

Institute for Safe Medication Practices (ISMP). Acute Care ISMP Medication Safety Alert! Ongoing risk: Transition adapters for ENFit syringes can defeat the purpose of ENFit itself. September 2017; Volume 22, Issue 19.

Pharmacy Syringes

- Metric scale gradients only
(Complies with ISMP Best Practices #5)
- Available in sizes 0.5 mL to 60 mL
- FDA 510(k) cleared for oral and enteral administration
- Sterile amber syringes available in sizes 1 to 12 mL



Self-righting, plugged tip caps promote aseptic technique and hands-free operation



Distinct gradient markings and bold numbers provide visual clarity

Amber Pharmacy Syringes

Legacy (Oral Tip)



Size	Sterile	Item #	Bulk	Item #
0.5 mL			BA-S05EO	16069
1 mL	BA-S1EO-S	16170	BA-S1EO	16070
3 mL	BA-S3EO-S	16171	BA-S3EO	16071
6 mL	BA-S6EO-S	16172	BA-S6EO	16072
12 mL	BA-S12EO-S	16173	BA-S12EO	16073
20 mL			BA-S20EO	16074
35 mL			BA-S35EO	16075
60 mL			BA-S60EO	16076

ENFit



Size	Sterile	Item #	Bulk	Item #
0.5 mL			BA-S05NC	816069
1 mL	BA-S1NC-S	816170	BA-S1NC	816070
3 mL	BA-S3NC-S	816171	BA-S3NC	816071
6 mL	BA-S6NC-S	816172	BA-S6NC	816072
12 mL	BA-S12NC-S	816173	BA-S12NC	816073
20 mL			BA-S20NC	816074
35 mL			BA-S35NC	816075
60 mL			BA-S60NC	816076

LDT

Blue Pharmacy Syringes



Legacy

Size	Bulk Product Code	Item #
0.5 mL	BB-S05EO	18077
1 mL	BB-S1EO	18078
3 mL	BB-S3EO	18079
6 mL	BB-S6EO	18080
12 mL	BB-S12EO	18081
20 mL	BB-S20EO	18082
35 mL	BB-S35EO	18083
60 mL	BB-S60EO	18084

Clear Pharmacy Syringes



Legacy

Size	Bulk Product Code	Item #
0.5 mL	BC-S05EO	14061
1 mL	BC-S1EO	14062
3 mL	BC-S3EO	14063
6 mL	BC-S6EO	14064
12 mL	BC-S12EO	14065
20 mL	BC-S20EO	14066
35 mL	BC-S35EO	14067
60 mL	BC-S60EO	14068

Purple Pharmacy Syringes



ENFit

Size	Bulk Product Code	Item #
0.5 mL	BP-S05NC	818077
1 mL	BP-S1NC	818078
3 mL	BP-S3NC	818079
6 mL	BP-S6NC	818080
12 mL	BP-S12NC	818081
20 mL	BP-S20NC	818082
35 mL	BP-S35NC	818083
60 mL	BP-S60NC	818084

LDT

Tamper-Evident Caps

A mechanical solution to help protect against diversion of controlled substances



Drug diversion concerns are driving the need to secure and protect syringe contents. NeoMed's Tamper-Evident Caps **enhance the security of the syringe and its contents** from the time it leaves the pharmacy until the time it is dispensed. Standardize your pharmacy with NeoMed Tamper-Evident Caps to ensure compatibility and protocol consistency.



Makes tampering immediately evident when the colored outer sleeve is broken



Tamper evidence protocols may help your pharmacy comply with DEA and USP guidelines



Available for both NeoMed Oral Tip Syringes and NeoConnect Syringes for ENFit

Legacy



Product Code	Item #
NM-5775	NM-5775

ENFit



Product Code	Item #
PNS-TENC	814775

Pharmacy Bottle Caps


- For medication filling into an ENFit syringe
- Designed for rapid filling operations
- Tethered cap for secure closure
- Designed to promote aseptic technique



ENFit



Size	Cap	Product Code	Item #
18 mm	A	PC-ANC	851001
20 mm	B	PC-BNC	851002
22 mm	C	PC-CNC	851003
24 mm	D	PC-DNC	851004
28 mm	E	PC-ENC	851005
33 mm	G	PC-GNC	851007
38 mm	I	PC-INC	851011



“We urge product vendors to make samples of ENFit syringes, caps, bottle adapters, feeding tubes, straws, and administration sets available as soon as possible.”

Institute for Safe Medication Practices Medication Safety Alert. ENFit Enteral Devices Are On Their Way...Important Safety Considerations for Hospitals. April 2015; 20;7.



DOSEMATE®



The **DoseMate** is an administration tip that is ergonomically designed for patient comfort for oral dosing directly to the mouth with an ENFit syringe



The **DoseMate DL** is ideal for dosing directly to the cheek or back of the mouth

ENFit



Product Code

Item #

POA-TNC

810105

ENFit



Product Code

Item #

POA-TNCL

810104

"...the syringe user should remove fluid that lies outside of the fluid path. The area between the male lumen and the outer ring (the "moat") is not part of the fluid path and should be free of fluid."

Global Enteral Device Supplier Association. GEDSA Position Statement Supporting ENFit® Low Dose Tip Syringe. August 2016.



Syringe-to-Syringe Coupler

Connects two syringes for transfer or mixing



ENFit



Purple	Item #	Product Code
Sterile Tray	810101	PCO-FNC
Bulk	810102	PCO-FNC-NS

Transfer Straws

For drawing up from a container into a syringe



Minimizes bacteria exposure when drawing formula or HBM from bottles



Can be used at bedside



Legacy



Orange	Item #	Product Code
6" Large Bore Straw	50038	NM-6ENSTB
Purple	Item #	Product Code
3" Small Bore Straw	51042	PNM-3ENST
6" Large Bore Straw	51038	PNM-6ENSTB

ENFit



Orange	Item #	Product Code
3" Small Bore Straw	850042	EXT-3STNC
6" Large Bore Straw	850038	EXT-6STBNC
Purple	Item #	Product Code
3" Small Bore Straw	851042	PEXT-3STNC
6" Large Bore Straw	851038	PEXT-6STBNC

ENFit Filtered Blunt Needle



ENFit



Item#	Length	Size
13061038	38 mm	18 ga

- Connects to an ENFit syringe to draw up medication from vials, glass ampules, and bottles
- Features a 5 µm filter for filling from glass ampules
- 18 gauge stainless steel needle with an angled tip allows for easy penetration through a vial's rubber membrane
- 38 mm long needle for easy use with most vials



ENFit



Item#	Product Code
851010	PC-EONC

Fill Cap Coupler

Allows ENFit syringes to connect with current oral press-in bottle adapters and other legacy bottle adapters

- Works with current oral bottle adapters: blue caps, press-in adapters and bung style adapters
- Promotes aseptic technique
- Easy to use for acute and home care settings
- Allows continued use of standard childproof bottle caps



NEOCONNECT[®] *at home*

Over-the-Counter Syringes

The NeoConnect[®] at home™ syringe line is ideal for patients in non-clinical and homecare settings.



Available for distribution or purchase without a prescription



Available for purchase through retail and online stores, DME's, and homecare companies



Size	Product Code	Item #
3 mL	Neo3	812002
6 mL	Neo6	812003
12 mL	Neo12	812004
60 mL	Neo60	812007



Syringes may be hand washed and reused up to **20 times** when cleaned per the cleaning instructions, allowing for a cost-effective and eco-friendly alternative to single use ENFit syringes.

NeoConnect at Home Syringes only available in the U.S.

® ENFit is a Registered Trademark or Trademark of Global Enteral Device Supplier Association, Inc. Used with permission. NeoConnect is a Registered Trademark of NeoMed, Inc. NeoMed is a Registered Trademark. Copyright 2020. NeoMed, Inc. All Rights Reserved. NM-1202 Rev 0.

 www.neomedinc.com 770.516.2225 info@neomedinc.com